

# High-Throughput Raman Process Monitoring of Downstream Protein Purification

Shaun Fraser, Shamus Driver, and Mark S. Kemper

Bruker Optics

40 Manning Road, Manning Park, Billerica, MA 01821, USA

## Raman Spectroscopy for Real-Time Bioprocess Monitoring

Raman spectroscopy is emerging as a powerful tool for in-line monitoring and downstream protein analysis. Enabled by HTVS™ technology, modern Raman systems deliver up to 30x higher signal, improving sensitivity, speed, and accuracy—even at low laser power. This allows real-time quantification of specific proteins and impurities during chromatographic separation, outperforming traditional UV-vis methods.

By providing actionable insights for process control, Raman supports PAT implementation, enhances product quality, and reduces downstream costs—unlocking significant gains in efficiency and profitability for biopharmaceutical manufacturing.

## Methods

An ÄKTA™ Start system was integrated with a Bruker HyperFlux™ Pro Plus Raman spectrometer using a flow cell and immersion probe for in-line measurements. Raman spectra were collected every 45 seconds (15 exposures, 3000 ms each), while UV absorbance at 260 nm correlated elution times.

### Measurements:

Raman spectra every 45 sec; UV absorbance at 260 nm for correlation.

### Chromatography:

HiPrep™ 16/60 Sephacryl™ S-200 HR column; buffer (25 mM sodium phosphate, 150 mM NaCl, pH 7.5).

### Samples:

25 mg BSA + 12.5 mg Cytochrome C in 5 mL buffer; injected at 0.5 mL/min; full elution after 95 mL




Fig. 2 (a) Raw Raman spectra of 5 mg/mL BSA in gel buffer (red), 2.5 mg/mL CytC in gel buffer (purple) and gel buffer (green). The region of interest (900-1500 cm<sup>-1</sup>) is further pre-processed and expanded in the inset (b), showing peak increases and shape changes in bands (amide II band ~1500 cm<sup>-1</sup>, amide III ~1300 cm<sup>-1</sup>) that are consistent with an increase of protein[6].

## Results

- Calibration Accuracy:** PLS model predicted BSA and CytC concentrations with errors of 0.109 mg/mL and 0.070 mg/mL (Fig. 2a & Fig. 2b).
- Elution Tracking:** Raman identified BSA peak at ~50 mL and CytC at ~82 mL, matching UV chromatogram (Fig. 3).
- Recovery:** Predicted totals closely matched injected amounts—93% for BSA and 89% for CytC (Fig. 4).
- Specificity:** Near-zero cross-contamination confirms high selectivity for individual proteins.

## Summary

Demonstrate HTVS™-enabled Raman spectroscopy for real-time, in-line quantification of individual proteins during downstream purification, compared to UV-vis

ÄKTA™ Start system integrated with Bruker HyperFlux™ Pro Plus Raman; monitored BSA and Cytochrome C elution using gel filtration and PLS calibration models (Fig. 1: Raman spectra; Fig. 2a & 2b: Calibration correlation).

Raman achieved high accuracy (errors: 0.109 mg/mL for BSA, 0.070 mg/mL for CytC), matched UV chromatogram, predicted recoveries of 93% and 89%, and showed strong specificity (Fig. 3: Elution spectra; Fig. 5: Predicted concentration profiles).

Enables faster, precise process control, optimizes protein collection, improves yield and purity, and supports PAT implementation for cost reduction and profitability

## Conclusion

- Raman predicts individual protein concentrations during purification.
- HTVS™ sensitivity enables fast, accurate measurements with superior SNR.
- Process Control: Optimizes collection, detects excursions, and improves yield.

Process Raman




Fig. 2 (a) Correlation plot resulting from the PLS model for the prediction of Cytochrome C constructed using the calibration samples. (b) Correlation plot resulting from the PLS model for the prediction of BSA.




Fig. 3 (Left) Rubber band subtraction, 19-pt smoothing and SNV corrected spectra of the gel filtration process colored by progressive cumulative flow volume.




Fig. 4 (Right) Predicted concentration chromatogram of BSA and CytC by volume. On the secondary axis are the units for the chromatogram generated by UV detection.